



The FHIR® Roadmap

ART-DECOR®

Tooling

dr Kai Heitmann, MD, FHL7

ART-DECOR® Expert Group
ADOT (ART-DECOR Open Tools)

Ministry of Health, Germany

info@kheitmann.de

Alexander Henket

ART-DECOR® Expert Group


Nictiz

HL7 The Netherlands


henket@nictiz.nl

ART-DECOR® Tool Suite...

- ...is an open-source tool and a methodology for various **multidisciplinary** stakeholders of healthcare information exchange
- ...supports comprehensive **collaboration of team members** within/between governance groups
- ...supports **creation, maintenance and federated sharing of HL7 templates** (HL7 STU), **value sets**, and more
- ...has a built-in **schematron engine** to allow easy **comprehensive validation and testing**




Supporting the “Life” Cycle






Functionality Overview

ART-DECOR® Project, Governance Group





Template Viewer / Editor

- Create / Maintain / Share  Templates + ...
- Template Viewer for HL7's STU format
- Two Template Editors
- Derive Artefacts: Documentation, Implementation Support, Validation/Testing

Demo 5: Vital Signs CDA Document Templates

DECOR Project Datasets Scenarios Terminology Templates Issues

Search

- CDA document level template
- Vital Signs Summary CDA document
- CDA header level template
 - CDA author
 - CDA typeId
 - CDA title
 - CDA confidentialityCode
 - CDA Organization Elements
 - CDA recordTarget
 - CDA id
 - CDA effectiveTime
 - CDA Person Elements
 - AssignedEntityElements
- CDA section level template
 - Vital Signs Section
 - CDA entry level template
 - Vital Sign Observation
 - Mean blood pressure
 - Body Weight
 - HearRatePeripheral

CDA recordTarget

Version: 2014-07-08 Status: Active
Version Label: Id: demo5-template-2001
Name: CDA recordTarget Display Name: CDA recordTarget

Description, model, context and examples

Description: The patient / client
Classification: cdaheaderlevel
Open/Closed: open

Used by / Uses

Used by Template Id	as	Name	as of
2.16.840.1.113883.3.1937.99.60.5.10.1	Inclusion	VitalSignsCDA document Vital Signs Summary CDA document	2014-07-08
Uses Template Id	as	Name	as of
2.16.840.1.113883.3.1937.99.60.5.10.900207	Inclusion	NationalPatientIdentifier Burgerservicenummer	2014-07-08

Relations







Type	Model / Template	as of
Specialization	2.16.840.1.113883.10.12.101	2005-09-07

Example


A simple patient, named John Doedidoe

```
<recordTarget>
  <patientRole>
    <id root="2.16.840.1.113883.2.4.6.3" extension="100202020"/>
    <name>
      <given>John</given>
      <family>Doedidoe</family>
    </name>
    <administrativeGenderCode code="M" codeSystem="2.16.840.1.113883.5.1"/>
    <birthTime value="19620219"/>
  </patientRole>
</recordTarget>
```

ART-DECORs FHIR® Roadmap


- FHIR® Profile/Resource **Viewer** 
- Multiple Exports from / Import to genuine ART-DECOR artefacts 
- Collaboration with FHIR® **Repositories**
 - simplifier.net 
- FHIR® **Terminology Services** 
 - Participation at the first Terminology Connectathon
- FHIR® **Implementation Guide Resources**
 - Using it for FHIR® and CDA IGs 
- FHIR® Profile **Editor**  not planned

ART-DECOR's FHIR® Roadmap



Logical <> StructureDefinition

Name	Flags	Card. Type	Description & Constraints
<div> <div>AllergyIntolerance</div> <div> <div>identifier</div> <div>onset</div> <div>recordedDate</div> <div>recorder</div> <div>patient</div> <div>reporter (Informant)</div> <div>substance</div> <div>status</div> <div>status-code-specification</div> </div> </div>		AllergyIntolerance	AllergyIntolerance Identifier StartDateTime RecordedDateTime Recorder Patient Informant CausativeSubstance Binding: VeroorzakendeStofCodeLijst (extensible) AllergyStatus Provides the specific CBB code for AllergyStatus URL: http://nictiz.nl/fhir/StructureDefinition/code-specification




Reference(nl-core-patient)
 Reference(nl-core-patient |
 nl-core-relatedperson | nl-
 core-practitioner)

0..1 Extension(CodeableConcept)

```

<element id="AllergyIntolerance.substance">
  <path value="AllergyIntolerance.substance"/>
  <short value="CausativeSubstance"/>
  <alias value="VeroorzakendeStof"/>
  <binding>
    <strength value="extensible"/>
    <valueSetReference>
      <reference value="http://decor.nictiz.nl/fhir/ValueSet/">
        <display value="VeroorzakendeStofCodeLijst"/>
      </valueSetReference>
    </binding>
    <mapping>
      <identity value="zib-allergieintolerantie"/>
      <map value="NL-CM-8.2.2"/>
    </mapping>
  </element>
  
```



Thus: things on the road

- More FHIR integration
- Terminology Services
- Issue and version management
- Diff functionality for ART-DECOR artifacts

- Organizational
 - ADOT: ART-DECOR Open Tools
- Question: FHIR validation and testing

Thank you

dr Kai Heitmann, MD, FHL7
ART-DECOR® Expert Group
ADOT (ART-DECOR Open Tools)
Ministry of Health, Germany
info@kheitmann.de

Alexander Henket
ART-DECOR® Expert Group
Nictiz
HL7 The Netherlands
henket@nictiz.nl

